

POSTER PRESENTATION

3D MR neurography targeted peripheral nerve ablation with MR-guided high intensity focused ultrasound (MR-HIFU): initial results of a feasibility study in a swine model

Robert Staruch^{1*}, Merel Huisman², Michelle Ladouceur-Wodzak³, Avneesh Chhabra³, Rajiv Chopra³

From Current and Future Applications of Focused Ultrasound 2014. 4th International Symposium Washington, D.C, USA. 12-16 October 2014

Background/introduction

MR-guided HIFU is an effective treatment for metastatic bone pain through periosteal nerve ablation, [1] and is being investigated for treating back pain through facet joint denervation.[2] For peripheral neuropathy, ultrasound-guided HIFU has been investigated preclinically as a means of achieving either an irreversible conduction block to treat severe spasticity [3], or a reversible partial conduction block to alleviate chronic pain.[4] However, ultrasound offers limited visualization of deeply situated pelvic nerves [5] and lacks the ability to measure thermal dose, which predicts the extent of changes in peripheral nerve histology and function [6]. Recently developed diffusion-prepared 3D MR neurography imaging techniques with fat suppression and nerve-selective T2-weighting [7] could improve targeting accuracy over ultrasound guidance. We present initial investigations into the use of MRI to guide HIFU ablation of peripheral nerves in a swine model. The objectives were 1) to evaluate the feasibility of identifying peripheral nerves using MR neurography on the clinical MR-HIFU system, 2) to monitor HIFU ablation of peripheral nerves using MR thermometry, and 3) to evaluate the ability to measure thermal lesions in peripheral nerves using contrast-enhanced T1-weighted images and thermal dose maps calculated from MR thermometry.

Methods

Experiments were approved by the local Institutional Animal Care and Use Committee. Volumetric MR-HIFU

¹Philips Research, Dallas, Texas, United States

was used to induce seven thermal lesions in the sciatic nerves of three pigs. 3D MR neurography and T1-weighted images at 3T were used for target identification and treatment planning. A single 8 or 12 mm treatment cell was used to cover the full width of each targeted nerve. Ultrasound exposures were performed under MR thermometry guidance in five image planes across the HIFU beam, and one plane along the beam axis. Sonications were performed at 1.2 MHz with acoustic power ranging from 160 to 300 W over fixed durations of 20 or 36 seconds (energy 3.2 to 10.8 kJ). Ablation dimensions were measured and compared using thermal dose maps, contrast-enhanced T1-weighted images, and gross pathology.

Results and conclusions

All targeted sciatic nerves were identifiable on MR neurography and T1-weighted images (Fig 1). For sonications at 160 to 300 W, MR thermometry measured peak temperatures of 60.3 to 85.7°C, with 240 equivalent minute thermal dose diameters of 8.5 to 15.9 mm (Fig 2). Thermal lesions were visible on late phase contrast-enhanced T1 (Fig 3), with dimensions matching the coagulated region observed at necropsy (Fig 4). Our preliminary results indicate that targeted peripheral nerve ablation is feasible with MR-HIFU. Diffusion-prep 3D MR neurography has potential for guiding therapy procedures where either nerve targeting or avoidance is desired.


Acknowledgements (Funding)


Funded provided by the Cancer Prevention and Research Institute of Texas and the M.R. and Evelyn Hudson Foundation. Robert Staruch is a paid employee of Philips Research. The authors thank Dave Hitt of Philips Healthcare for technical support in customizing the MR neurography protocol for the swine model.




© 2015 Staruch et al; licensee BioMed Central Ltd. This is an Open Access article distributed under the terms of the Creative Commons Attribution License (http://creativecommons.org/licenses/by/4.0), which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited. The Creative Commons Public Domain Dedication waiver (http:// creativecommons.org/publicdomain/zero/1.0/) applies to the data made available in this article, unless otherwise stated.


Full list of author information is available at the end of the article














Authors' details

¹Philips Research, Dallas, Texas, United States. ²University Medical Center Utrecht, Utrecht, Netherlands. ³University of Texas Southwestern Medical Center, Dallas, Texas, United States.

Published: 30 June 2015

References

- Hurwitz , et al: Magnetic Resonance-Guided Focused Ultrasound for Patients With Painful Bone Metastases: Phase III Trial Results. JNCI 2014, 106(5).
- Weeks, et al: MRI-guided focused ultrasound (MRgFUS) to treat facet joint osteoarthritis low back pain-case series of an innovative new technique. Eur Radiol 2012, 22:2822-35.
- 3. Foley, *et al*: Image-guided high-intensity focused ultrasound for conduction block of peripheral nerves. *Ann Biomed Eng* 2007, **35**:109-19.
- 4. Foley, *et al*: Effects of high-intensity focused ultrasound on nerve conduction. *Muscle Nerve* 2008, **37**:241-50.
- Fritz, et al: Magnetic resonance neurography-guided nerve blocks for the diagnosis and treatment of chronic pelvic pain syndrome. *Neuroimag Clin N Am* 2014, 24(1):211-34.
- Vujaskovic, et al: Effects of intraoperative hyperthermia on peripheral nerves: neurological and electrophysiological studies. Int J Hyperthermia 1994, 10(1):41-9.
- 7. Yoneyama, et al: Rapid high resolution MR neurography with a diffusionweighted pre-pulse. Magn Reson Med Sci 2013, 12(2):111-9.

doi:10.1186/2050-5736-3-S1-P74

Cite this article as: Staruch *et al.*: **3D MR neurography targeted** peripheral nerve ablation with MR-guided high intensity focused ultrasound (MR-HIFU): initial results of a feasibility study in a swine model. *Journal of Therapeutic Ultrasound* 2015 **3**(Suppl 1):P74.

Submit your next manuscript to BioMed Central and take full advantage of:

- Convenient online submission
- Thorough peer review
- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar
- Research which is freely available for redistribution

Submit your manuscript at www.biomedcentral.com/submit

BioMed Central