

Open Access

POSTER PRESENTATION

Thermal ablation of a confluent lesion in the porcine kidney with magnetic resonance guided high intensity focused ultrasound

Johanna van Breugel^{1*}, Joost Wijlemans¹, Martijn de Greef¹, Gerald Schubert², Maurice van den Bosch¹, Chrit Moonen¹, Mario Ries¹

From Current and Future Applications of Focused Ultrasound 2014. 4th International Symposium Washington, D.C, USA. 12-16 October 2014

Background/introduction

Since approximately 1.6 percent of men and women will be diagnosed with kidney and renal pelvis cancer during their lifetime, there is a growing interest in non-invasive kidney sparing therapy for renal cancer. As a consequence, several patient studies investigated the feasibility of high intensity focused ultrasound (HIFU) for the thermal ablation of renal masses. The majority of these studies used either a hand-held extracorporeal ultrasound transducer with ultrasound imaging for guidance or a laparoscopic approach. Drawbacks of these techniques are the lack of respiratory motion compensation, no means to observe the energy deposition in real time, the complexity of the probe positioning, and the risks of bleeding and tumor spillage. Alternatively, recent preclinical studies have demonstrated the feasibility of magnetic resonance guided high intensity focused ultrasound (MR-HIFU) interventions on the kidney with respect to motion compensated real-time thermometry and acoustic energy delivery. Here, we extend this prior work to investigate in an animal study if MR-HIFU can deliver a reliable confluent volumetric lesion in the renal cortex in a clinically relevant time-frame.

Methods

An anesthetized Dalland land pig was placed on its right side on a clinical Sonalleve MR-HIFU therapy system, which is integrated with a 1.5T Achieva MRI (Philips Healthcare) with minor modifications. Both acoustic energy delivery and MR-thermometry were respiratory gated (3mm gating window, ~70% duty cycle) and active

¹University Medical Center Utrecht, Utrecht, Netherlands Full list of author information is available at the end of the article




Figure 1 Planning of 7 treatment cells in the cortex of the kidney

surface cooling was employed to prevent undesired nearfield damage. A honeycomb pattern of seven ablation cells (12-17s, 450W acoustic power, 4x4x10 mm3) was positioned in the cortex of the kidney as shown in figure 1. The therapeutic endpoint was evaluated non-invasively at the end of the intervention. Hereto, lethal thermal dose



© 2015 van Breugel et al; licensee BioMed Central Ltd. This is an Open Access article distributed under the terms of the Creative Commons Attribution License (http://creativecommons.org/licenses/by/4.0), which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited. The Creative Commons Public Domain Dedication waiver (http://creativecommons.org/publicdomain/zero/1.0/) applies to the data made available in this article, unless otherwise stated.

estimates based on MR-thermometry and a non-perfused volume (NPV) measurement using dynamic contrast enhanced T1-weighted MRI (DCE-MRI) were performed. Subsequently, the animal was euthanized and the extent of the induced necrosis examined using a cellular viability staining (nicotinamide adenine dinucleotide, NADH).

Figure 2 Contrast enhanced MR scan shows a non-perfused

Results and conclusions

volume in the kidney

All ablation cells reached peak temperatures > 80°C during the sonications, resulting in a lethal thermal dose over the entire delineated target area. DCE-MRI displayed a confluent non-perfused volume within the cortex and partly within the medulla with a volume of ± 2 ml as shown in figure 2. NADH staining reconfirmed a confluent nonviable volume of approximately (12 x 16 x 10 mm3) (Figure 3). No undesired tissue damage in adjacent areas has been observed. These first results indicate that current MR-guided clinical HIFU equipment might be suitable for non invasive therapy of renal masses. Future work will

need to demonstrate the reproducibility of the findings and investigate potential adverse effects (cutaneous and subcutaneous damage) as a preparation for a clinical study.

Acknowledgements (Funding)

This study was performed within the framework of CTMM, the Center for Translational Molecular Medicine, project VOLTA (grant 05T-201).

Authors' details

¹University Medical Center Utrecht, Utrecht, Netherlands. ²Philips Healthcare, Vantaa, Finland.

Published: 30 June 2015

doi:10.1186/2050-5736-3-S1-P76 **Cite this article as:** van Breugel *et al.*: **Thermal ablation of a confluent** lesion in the porcine kidney with magnetic resonance guided high intensity focused ultrasound. Journal of Therapeutic Ultrasound 2015 3 (Suppl 1):P76.

van Breugel et al. Journal of Therapeutic Ultrasound 2015, 3(Suppl 1):P76

http://www.jtultrasound.com/content/3/S1/P76

Figure 3 NADH staining shows a non-viable area (yellowish) within viable kidney tissue (blue).

